

## T-Ring™ Tourniquet for Digits

PLS1810 Bleeding Control, Maximum Application Time of up to 90 Minutes, 25 per box

Stop the blood flow from a deep finger wound with the exact-right amount of pressure using this physician-designed finger tourniquet ring.

- Stops the bleeding immediately
- Applies the ideal safe amount of pressure, every time, especially compared to other methods where excess pressure risks nerve damage
- Designed by a physician
- Two-point discriminator assesses nerve function
- This one-size finger-and-toe tourniquet fits and works correctly on all digit sizes, from pinky to big toe
- Simple slide-on application; no tying, cutting, rolling, clamping or guessing
- Easily "stretches" to apply over skin flaps or deformities
- Break-away design makes removal easy and painless compared to sliding back off
- Highly visible to prevent forgetting that it's in place
- Each unit is latex free, individually packaged and sterilized

### Specifications

Style	Slip-On
Dimensions	1.438" W x 1.875" L
Application	Bleeding Control
Brand	T-Ring
Fits	One Size Fits All
Intended Use	Maximum Application Time of up to 60 Minutes
Color	Red/Silver
Special Feature	Breakaway Exterior; Attached Discriminator
Sold as	25 per box
Weight	.37 lbs.
Composition	Polycarbonate Thermoplastic elastomers
Shelf Life	5 years
UNSPSC	41104104

## Metric Equivalent

Weight	168 g
Dimensions	3.7cm W x 4.8cm L

## Technical Information



**One Pork Avenue • Tipton, PA 16684-0304**  
1-855-493-4647 • Fax: 1-800-621-7447 • [newpig.com](http://newpig.com) • [hothogs@newpig.com](mailto:hothogs@newpig.com)